



For a Better Life

# IN-VITRO DIAGNOSTICS PRODUCT RANGE

**2024** Edition I

[www.ameworldwide.com](http://www.ameworldwide.com)  
[www.almaniads.com](http://www.almaniads.com)

## Catalog


New Product	1-2
Fertility&Eugenic Test	3-4
Infectious Disease Test	5-6
Drug Of Abuse Test	7-8
SARS-CoV-2 Rapid Test Kit	9
Product Overview	10-13



For a Better Life

www.ameworldwide.com  
www.almaniads.com

**New Product**



**Urine Sample HPV Virus Staining Kit**

**H.Pylori Antigen Test By Dental Deposit**

**Features:**

- ✓ **Specimen type: Urine/Dental Deposit**
- ✓ **Self-testing at home**
- ✓ **Testing time: 15 minutes**



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# Fertility Products

Product Name	Format	Specimen	Sensitivity	Reg.	
Pregnancy(HCG) Test	Strip/Cassette/Midstream	Urine	25/20/10mIU/ml	K123050* K232864	CE <sub>0123</sub> FSC Anvisa EFDA
Pregnancy(HCG) Test	Strip/Cassette	Urine/Serum	25mIU/ml		FSC
Ovulation(LH) Test	Strip/Cassette/Midstream	Urine	25mIU/ml	D163970	CE <sub>0123</sub> FSC
FSH Test	Strip/Cassette/Midstream	Urine	25mIU/ml		CE <sub>0123</sub> FSC
SPERM Concentration Test	Cassette	Seminal fluid	15x10 <sup>6</sup> /ml		CE <sub>0123</sub> FSC
Digital Pregnancy Test	Midstream	Urine	25mIU/ml		CE <sub>0123</sub> FSC

\*K123050 only sensitivity 25mIU/ml approved



## Electrical midstream for gestational weeks



Show the number of weeks pregnant  
Double confirmation



Accurate results from the estimated period of menstruation  
Get the results four days in advance



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# Infectious Disease




Dengue IgM/IgG



Dengue NS1



Dengue Combo



### Dengue IgM and/or IgG Rapid Test

Format:Cassette

Specimen:Serum(Plasma)/Whole Blood

Dengue IgM and/or IgG Rapid Test Cassette is a lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of IgM and IgG anti-dengue in human serum (or plasma). The kit is a diagnostic device for qualitative detection .



Pf/pv



Pf



Pan/Pf



HIV1+2



Negative:



Positive:

P. falciparum malaria



Invalid:



P.vivax malaria:



P.vivax malaria and P. falciparum malaria:



TP



HCV

Format:Cassette

Specimen:Whole Blood

For the rapid qualitative determination of Malaria P.falciparum specific histidine rich protein-2(Pf HRP-2)and malaria P.vivax specific lactate dehydrogenase (pLDH) in human blood as an aid in the diagnosis of Malaria infection.



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# INFECTIOUS DISEASE

## Tropics Infectious Disease Rapid Test Kits

Product Name	Format	Specimen	Reg.
Malaria Pf	Cassette	WB	CFDA FSC EFDA
Malaria Pv+Pf	Cassette	WB	CFDA FSC EFDA
Malaria Pan+Pf	Cassette	WB	FSC EFDA
Dengue IgG/IgM	Cassette	S/P/WB	FSC
Dengue NS1 antigen	Cassette	S/P/WB	FSC
Dengue IgG/IgM NS1 combo	Cassette	S/P/WB	FSC
Typhoid IgG/IgM	Cassette	S/P/WB	FSC

## Infectious Disease Rapid Test Kits

Product Name	Format	Specimen	China Reg.No
HAV IgM	Cassette	S/P/WB	CFDA FSC
HEV IgM	Cassette	S/P/WB	CFDA FSC
HAV IGM HEV IGM COMBO	Cassette	S/P/WB	CFDA FSC
HBsAg	Strip/Cassette	S/P/WB	CFDA FSC EFDA
HBsAb	Strip/Cassette	S/P/WB	CFDA FSC EFDA
HbeAg	Strip/Cassette	S/P/WB	FSC EFDA
HBeAb	Strip/Cassette	S/P/WB	FSC EFDA
HBcAb	Strip/Cassette	S/P/WB	FSC EFDA
HBcAb IgM	Strip/Cassette	S/P/WB	FSC EFDA
HBV COMBO	Cassette	S/P/WB	FSC EFDA
HCV	Strip/Cassette	S/P/WB	CFDA FSC EFDA
TB	Cassette	S/P/WB	FSC
SYPHILIS	Strip/Cassette	S/P/WB	CFDA FSC EFDA
HIV 1+2	Strip/Cassette	S/P/WB	CFDA FSC
HIV 1+2 TRILINES	Cassette	S/P/WB	FSC
HIV 1+2+O	Cassette	S/P/WB	FSC
TOX-IgM	Cassette	S/P/WB	FSC
RV-IgM	Cassette	S/P/WB	FSC
CMV-IgM	Cassette	S/P/WB	FSC
HSV- I -IgM	Cassette	S/P/WB	FSC
HSV- II -IgM	Cassette	S/P/WB	FSC
TOX-IgG	Cassette	S/P/WB	FSC
RV-IgG	Cassette	S/P/WB	FSC
CMV-IgG	Cassette	S/P/WB	FSC
HSV- I -IgG	Cassette	S/P/WB	FSC
HSV- II -IgG	Cassette	S/P/WB	FSC
TORCH-IgM(COMBO-5)	Cassette	S/P/WB	FSC
TORCH-IgG(COMBO-5)	Cassette	S/P/WB	FSC
H.pylori Antigen	Strip/Cassette	Feces	FSC EFDA
H.pylori Antibody	Strip/Cassette	S/P/WB	CE CFDA FSC EFDA

\*CFDA: Registered in China

\*EFDA: Registered in Ethiopia





For a Better Life

# DRUG OF ABUSE KITS




## Features & Benefits

- Easy to use
- Results in 5 minutes
- Visually interpreted
- Multiple formats: dipcard, cassette, urine cup, saliva device
- High accuracy
- Cost effective

### Strip

- Immerse the strip into urine
- Read results at 5 minutes




### Cassette

- Add 2-3 full drops (approx. 100uL) of urine
- Read results at 5 minutes



### Panel

- Immerse the strip into urine
- Read results at 5 minutes



### Cup

- Collect the fresh urine sample in the cup
- Read results at 5 minutes





For a Better Life

www.ameworldwide.com  
www.almaniads.com

# DRUG OF ABUSE KITS

## STRIPS, PANELS AND CUPS

Product	Method	Specimen	Sensitivity	Certificate
Morphine (MOP, OPI)	Colloidal Gold	Urine	300ng/ml	K151557 FSC CE
Methyl amphetamine (MET)	Colloidal Gold	Urine	1000ng/ml	K151557 FSC CE
Amphetamine(AMP)	Colloidal Gold	Urine	1000ng/ml	FSC CE
Cocaine (COC)	Colloidal Gold	Urine	300ng/ml	FSC CE
Ketamine (KET)	Colloidal Gold	Urine	1000ng/ml	FSC CE
Methadone (MTD)	Colloidal Gold	Urine	300ng/ml	FSC CE
Marijuana(THC)	Colloidal Gold	Urine	50ng/ml	K152643 FSC CE
Barbital (BAR)	Colloidal Gold	Urine	300ng/ml	FSC CE
MDMA	Colloidal Gold	Urine	500ng/ml	K152643 FSC CE
Benzodiazepines (BZO)	Colloidal Gold	Urine	300ng/ml	FSC CE
Bupivacaine (BUP)	Colloidal Gold	Urine	10ng/ml	FSC CE
Phencyclidine (PCP)	Colloidal Gold	Urine	25 ng/ml	FSC CE
Oxycodone(OXY)	Colloidal Gold	Urine	100 ng/ml	FSC CE
TCA	Colloidal Gold	Urine	1000 ng/ml	FSC CE
Tramadol(TRA)	Colloidal Gold	Urine	100 ng/ml	FSC CE
EDDP	Colloidal Gold	Urine	100 ng/ml	FSC CE
Propoxyphen(PPX)	Colloidal Gold	Urine	300 ng/ml	FSC CE
Caffeine(CAF)	Colloidal Gold	Urine	8ug/ml	FSC CE
Methaqualone(MQL)	Colloidal Gold	Urine	300 ng/ml	FSC CE
Fentanyl(FYL)	Colloidal Gold	Urine	200 ng/ml	FSC CE
Synthetic Cannabis(K2)	Colloidal Gold	Urine	50 ng/ml	FSC CE
Carisoprodol/Somadril(SOMA)	Colloidal Gold	Urine	1000 ng/ml	FSC CE
Morphine (MOP, OPI)	Colloidal Gold	Saliva	15 ng/ml	FSC CE
Methyl amphetamine (MET)	Colloidal Gold	Saliva	50 ng/ml	FSC CE
Ketamine (KET)	Colloidal Gold	Saliva	100 ng/ml	FSC CE
Multi-Drug 2-12 Test Panel	Colloidal Gold	Urine	--	FSC CE
Multi-Drug 2-15 Test Cup	Colloidal Gold	Urine	--	FSC CE



Cassette



Strip



Panel








For a Better Life

www.ameworldwide.com  
www.almaniads.com


# SARS-CoV-2 TEST KIT

## SARS-CoV-2 Antigen Rapid Test


### SARS-CoV-2 Antigen Rapid Test



### SARS-CoV-2 Antigen Test for Professional Use



### Antigens Multiplex (FluA/FluB/SARS-CoV-2) Rapid Test



Product	Specimen	Packing	Reg.
SARS-CoV-2 Antigen	Oropharyngeal Nasopharyngeal	1pcs/Box,5pcs/Box,25pcs/Box	CE,Anvisa,Bfarm
SARS-CoV-2 Antigen(saliva)	Saliva	1pcs/Box,5pcs/Box,25pcs/Box	CE,Bfarm
FluA/FluB/SARS-Cov-2 Antigens Multiplex	Oropharyngeal Nasopharyngeal	1pcs/Box,25pcs/Box	CE,Anvisa
SARS-CoV-2 IGG/IgM	WB/S/P	1pcs/Box,25pcs/Box	CE,Anvisa
SARS-CoV-2 Multiplex Antibody	WB/S/P	25pcs/Box	CE
SARS-CoV-2 Neutralizing Antibody	WB/S/P	25pcs/Box	CE



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# PRODUCT OVERVIEW

## Fertility&Eugenic Test Kit

Product	Specimen	Format	Catalog
HCG Pregnancy	Urine	Strip/Cassette/Midstream	FE04-01/02/03
HCG Pregnancy	Serum/Urine	Strip/Cassette	FE04-01/02/04
LH ovulation	Urine	Strip/Cassette/Midstream	FE05-01/02/03
fetal fibronectin (fFN)	cervical secretions	Cassette	FE02-01
Follicle-Stimulating Hormone(FSH)	Urine	Strip/Cassette	FE06-01/02
TSH	WB/S/P	Strip/Cassette	FE0701/02
SPERM Concentration	Sperm	Cassette	FE08-01/02
Rv-IgM	WB/S/P	Strip/Cassette	FE09-01/02
CMV-IgM	WB/S/P	Strip/Cassette	FE10-01/02
TOXO IgM	WB/S/P	Strip/Cassette	FE11-01/02
HSV-I-IgM	WB/S/P	Strip/Cassette	FE12-01/02
HSV-II-IgM	WB/S/P	Strip/Cassette	FE12-03/04
Torch-IgM ( 5 in 1)	WB/S/P	Combo(5 in 1)	FE13-01

## Infectious Disease Test

Product	Specimen	Format	Catalog
SARS-CoV-2 Antigen	Oropharyngeal Nasopharyngeal	Cassette	COV01-01
SARS-CoV-2 Antigen(saliva)	Saliva	Cassette	COV02-01
FluA/FluB/SARS-Cov-2 Antigens Multiplex	Oropharyngeal Nasopharyngeal	Cassette	COV03-01
SARS-CoV-2 IGG/IgM Antibody	WB/S/P	Cassette	COV04-01
SARS-CoV-2 Multiplex Antibody	WB/S/P	Cassette	COV05-01
SARS-CoV-2 Neutralizing Antibody	WB/S/P	Cassette	COV06-01
Mycoplasma Pneumoniae(MP) IgM	WB/S/P	Cassette	ID01-01
InfluenzaA /InfluenzaB (FluA+B) Antigens	Oropharyngeal Nasopharyngeal	Cassette	ID02-01
Rotavirus And Adenovirus Antigen combo detection	Feces	Combo(2 in 1)	ID04-01
HAV-IgM	WB/S/P	Strip/Cassette	ID05-01/02
HAV/HEV IgM	WB/S/P	Combo(2 in 1)	ID06-01
HBsAg	WB/S/P	Strip/Cassette	ID07-01/02
HBV Combo(HBsAg,HBsAb,HBeAg,HBeAb,HBcAb)	WB/S/P	Combo(5 in 1)	ID08-01
HCV	WB/S/P	Strip/Cassette	ID09-01/02
HEV-IgM	WB/S/P	Strip/Cassette	ID10-01/02



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# PRODUCT OVERVIEW

HIV 1/2 Antibody	WB/S/P	Strip/Cassette	ID11-01/02
HIV 0/1/2 triline	WB/S/P	Cassette	ID11-03
Syphilis/TP	WB/S/P	Strip/Cassette	ID12-01/02
Malaria PF	WB	Cassette	ID13-01
Malaria PF/Pv	WB	Cassette	ID14-01
Malaria PF/Pan	WB	Cassette	ID15-01
Dengue IgG/IgM	WB/S/P	Cassette	ID16-01
Dengue NSI Antigen	WB/S/P	Cassette	ID17-01
Typhoid IgG/IgM	WB/S/P	Cassette	ID18-01
Chlamydia	Cervical Mucus	Cassette	ID19-01
Gonorrhea	Cervical Mucus	Cassette	ID20-01
H.pylori Antigen	Feces	Strip/Cassette	ID21-01
H.pylori Antibody	WB/S/P	Strip/Cassette	ID22-01
Urine Sample HPV Virus Staining Kit	Urine	Cassette	
H.Pylori Antigen Test By Dental Deposit	Dental Deposit	Cassette	

## Tumor marker Test

Product	Specimen	Format	Catalog
FOB(Fecal Occult Blood)	Feces	Strip/Cassette	TM01-01/02
CEA(Carcino-Embryonic Antigen)	WB/S/P	Strip/Cassette	TM02-01/02
AFP(Alpha Fetoprotein)	WB/S/P	Strip/Cassette	TM02-03/04
PSA(Prostate Specific Antigen)	WB/S/P	Strip/Cassette	TM02-04/04

## Cardiac Marker Test

Product	Specimen	Expiry	Catalog
Troponin(cTnI) Quantitative	Feces	24 Monthes	CM01-01
NT-proBNP Quantitative	WB/S/P	24 Monthes	CM02-01
Myohemoglobin(Myo) Quantitative	WB/S/P	24 Monthes	CM03-01
CK-MB Quantitative	WB/S/P	24 Monthes	CM04-01
H-FABP Quantitative	WB/S/P	24 Monthes	CM05-01
Troponin(cTnI) I	WB/S/P	18 Monthes	CM06-01
cTnI/Myo/CK-MB 3 in 1	WB/S/P	18 Monthes	CM07-01



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# PRODUCT OVERVIEW

## Drug Of Abuse Test

Product	Specimen	Format	Catalog
Morphine (MOP, OPI)	Urine	Strip/Cassette	DOA01-01/02
Methyl amphetamine (MET)	Urine	Strip/Cassette	DOA02-01/02
Amphetamine(AMP)	Urine	Strip/Cassette	DOA03-01/02
Cocaine (COC)	Urine	Strip/Cassette	DOA04-01/02
Ketamine (KET)	Urine	Strip/Cassette	DOA05-01/02
Methadone (MTD)	Urine	Strip/Cassette	DOA06-01/02
Marijuana(THC)	Urine	Strip/Cassette	DOA07-01/02
Barbital (BAR)	Urine	Strip/Cassette	DOA08-01/02
MDMA	Urine	Strip/Cassette	DOA09-01/02
Benzodiazepines (BZO)	Urine	Strip/Cassette	DOA10-01/02
Bupivacaine (BUP)	Urine	Strip/Cassette	DOA11-01/02
Phencyclidine (PCP)	Urine	Strip/Cassette	DOA12-01/02
Oxycodone(OXY)	Urine	Strip/Cassette	DOA13-01/02
TCA	Urine	Strip/Cassette	DOA14-01/02
Tramadol(TRA)	Urine	Strip/Cassette	DOA15-01/02
Synthetic Cannabis(K2)	Urine	Strip/Cassette	DOA16-01/02
Carisoprodol/Somadril(SOMA)	Urine	Strip/Cassette	DOA17-01/02
Propoxyphen(PPX)	Urine	Strip/Cassette	DOA18-01/02
Caffeine(CAF)	Urine	Strip/Cassette	DOA19-01/02
Methaqualone(MQL)	Urine	Strip/Cassette	DOA20-01/02
EDDP	Urine	Strip/Cassette	DOA21-01/02
Bupivacaine (BUP)	Urine	Strip/Cassette	DOA22-01/02
Multi-Drug 2-10 Test Panel	Urine	Combo	DOA23
Multi-Drug 2-15 Test Cup	Urine	Cup	DOA24
Morphine (MOP, OPI)	Saliva	Cassette	DOA25-01
Methyl amphetamine (MET)	Saliva	Cassette	DOA26-01
Ketamine (KET)	Saliva	Cassette	DOA27-01
MOR/MET test combo	Saliva	Cassette	DOA28
MOR/MET/KET test combo	Saliva	Cassette	DOA29



For a Better Life

www.ameworldwide.com  
www.almaniads.com

# PRODUCT OVERVIEW

## Pet Detection

Feline 🐾	Canine 🐾
Feline Pregnancy (FCG)	Canine Distemper Ag (CDV Ag)
Feline Leishmaniasis (FLsh)	Anaplasma IgG (ANA IgG)
Feline Panleukopenia (FPL)	Canine Infectious Hepatitis Ag (CAV Ag)
Feline Anaplasma IgG (FANA IgG)	Canine Heartworm AgTest (CHW Ag)
Feline Immuno deficiency IgG (FIV IgG)	CPV/CCV/Giardia Combo (CPV-CCV-Gia)
Feline Leukemia Ag (FELV Ag)	Canine Parvo/Corona Ag Combo (CPV-CCV Ag)
Feline Calicivirus Ag (FCV Ag)	CPV/CDV Ag Combo (CPV-CDV Ag)
Feline Infectious Peritonitis IgG (FIPV IgG)	Canine Parvo Ag (CPV Ag)
Feline Coronavirus Ag (FCoV Ag)	Canine Pregnancy Relaxin (RLN)
FelinePanleukopenia/Corona/Giardia Combo (FPV-FCOV-GIA)	Toxoplasma IgG/IgM IgG (Toxo IgG)
Feline Panleucopenia Ag (FPV Ag)	Canine Giardia Ag (Gia Ag)
FIV IgG/FeLV Ag Combo (FIV IgG-FELV Ag)	Brucella (Bru)
Canine 🐾	Canine 🐾
Rlgies IgG	Canine Distemper/Adeno/Influenza Combo (CDV-CAV-CIV)
Rlgies Ag (FPV Ag)	Canine Distemper/Adeno Ag Combo (CDV-CAV)
Ehrlichia/Lyme/Anaplasma/Heartorm Combo (EHR-LYM-ANA-CHW)	
Ehrlichia/Lyme/Anaplasma Combo (EHR-LYM-ANA)	
Ehrlichia/Anaplasma Combo (EHR-ANA)	
Ehrlichia canis IgG (E.canis IgG)	
Canine Brucella IgG (C.Bru IgG )	
Canine Coronavirus Ag (CCV Ag )	
Lyme Disease IgG (Lyme IgG)	
Canine Influenza Ag Rapid (CIV Ag)	
Leishmania canis IgG (LSH IgG)	





## CONTACT US



For a Better Life

**GLOBAL MARKETING AUTHORIZATION HOLDER**  
**ALMANIA DRUGS STORE**  
Sharjah, United Arab Emirates  
[www.almaniads.com](http://www.almaniads.com)  
[www.ameworldwide.com](http://www.ameworldwide.com)

North America Distribution Office:  
**AME WORLDWIDE L.L.C.**  
New Jersey, United States of America  
[www.ameworldwide.com](http://www.ameworldwide.com)

